

Product Name: STAT2 (5I18) Rabbit Monoclonal Antibody
Catalog #: AMRe18349

Summary

Production Name	STAT2 (5I18) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,FC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
Purification	Affinity purification

Immunogen

Gene Name	STAT2
Alternative Names	
Gene ID	6773.0
SwissProt ID	P52630.

Application

Dilution Ratio	WB 1:1000-1:5000, IHC-P/IF-P 1:100-1:200, ICC/IF 1:50, FCM 1:50-1:100
Molecular Weight	98 kDa


Product Name: STAT2 (5I18) Rabbit Monoclonal Antibody
Catalog #: AMRe18349




Background

Research Area


Image Data



Western blot analysis of STAT2 in K562 lysates using STAT2 antibody.




Immunohistochemistry analysis of paraffin-embedded Human lung cancer using STAT2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Product Name: STAT2 (5I18) Rabbit Monoclonal Antibody
Catalog #: AMRe18349



Immunofluorescence analysis of STAT2 in A431 using STAT2 antibody.



Western blot analysis of STAT2 in HeLa , Rat brain lysates using STAT2 antibody.

Note

For research use only.