
Product Name: STC1 (4H6) Rabbit Monoclonal Antibody**Catalog #: AMRe18377**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1000-1:5000
Molecular Weight	28kDa

Antigen Information

Gene Name	STC1
Alternative Names	Stanniocalcin; STC; STC-1; Stc1;
Gene ID	6781.0
SwissProt ID	P52823
Immunogen	A synthetic peptide of human Stanniocalcin 1


Background

Stimulates renal phosphate reabsorption, and could therefore prevent hypercalcemia. Stimulates renal phosphate

reabsorption, and could therefore prevent hypercalcemia.

Research Area

Image Data



Western blot analysis of STC1 expression in Human Prostate lysate.