
Product Name: Syntaxin 3 (18G14) Rabbit Monoclonal Antibody**Catalog #: AMRe18513**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1000-1:5000,IHC 1:50-1:100,ICC/IF 1:50-1:100
Molecular Weight	33kDa

Antigen Information

Gene Name	STX3
Alternative Names	Stx3; STX3A; Syntaxin 3A;
Gene ID	6809.0
SwissProt ID	Q13277
Immunogen	A synthetic peptide of human Syntaxin 3

Background


Potentially involved in docking of synaptic vesicles at presynaptic active zones. Potentially involved in docking of synaptic

vesicles at presynaptic active zones.

Research Area

Neuroscience

Image Data



Western blot analysis of Syntaxin 3 expression in HepG2 cell lysate.