

Summary

Production Name	TAB1 (1R7) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,FC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	TAB1
Alternative Names	MAP3K7IP1; TAB1; TAK1 binding protein 1;
Gene ID	10454.0
SwissProt ID	Q15750.

Application

Dilution Ratio	WB 1:1000, IHC-P/IF-P 1:100-1:200, ICC/IF 1:100-1:200, FCM 1:20
Molecular Weight	55kDa

Background


Product Name: TAB1 (1R7) Rabbit Monoclonal Antibody
Catalog #: AMRe18594



May be an important signaling intermediate between TGF β receptors and MAP3K7/TAK1. May play an important role in mammalian embryogenesis. May be an important signaling intermediate between TGF β receptors and MAP3K7/TAK1. May play an important role in mammalian embryogenesis.

Research Area

Image Data



Western blot analysis of TAB1 expression in MCF-7 cell lysate.

Note

For research use only.