

Product Name: TRAP (12F3) Rabbit Monoclonal Antibody
Catalog #: AMRe19221

Summary

Production Name	TRAP (12F3) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,IP,IF-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	ACP5
Alternative Names	TRAP; SPENCDI; TrATPase;
Gene ID	54.0
SwissProt ID	P13686.

Application

Dilution Ratio	WB 1:1000-1:5000, IHC-P/IF-P 1:100, ICC/IF 1:100, IP 1:20
Molecular Weight	37kDa


Product Name: TRAP (12F3) Rabbit Monoclonal Antibody
Catalog #: AMRe19221

Background

Involved in osteopontin/bone sialoprotein dephosphorylation. Its expression seems to increase in certain pathological states such as Gaucher and Hodgkin diseases, the hairy cell, the B-cell, and the T-cell leukemias. Involved in osteopontin/bone sialoprotein dephosphorylation. Its expression seems to increase in certain pathological states such as Gaucher and Hodgkin diseases, the hairy cell, the B-cell, and the T- cell leukemias.

Research Area

Image Data



Western blot analysis of Tartrate Resistant Acid Phosphatase expression in A431 cell lysate.

Note

For research use only.