

Product Name: JAK1 Rabbit Monoclonal Antibody
Catalog #: AMRe21094



Summary

Production Name	JAK1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC,IF,IP,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG,Kappa
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
Purification	Protein A

Immunogen

Gene Name	JAK1 JAK1A JAK1B
Alternative Names	Tyrosine-protein kinase JAK1;Janus kinase 1;JAK-1;
Gene ID	3716.0
SwissProt ID	P23458.

Application

Dilution Ratio	IHC 1:100-1:1000;WB 1:2000-1:10000;IF 1:200-1:1000;ELISA 1:5000-1:20000;IP 1:50-1:200;
Molecular Weight	Calculated MW:133kD;Observed MW:133kD

Background

Cell localization:Endomembrane system; Peripheral membrane protein. Wholly intracellular, possibly membrane


Product Name: JAK1 Rabbit Monoclonal Antibody
Catalog #: AMRe21094



associated..This gene encodes a membrane protein that is a member of a class of protein-tyrosine kinases (PTK) characterized by the presence of a second phosphotransferase-related domain immediately N-terminal to the PTK domain. The encoded kinase phosphorylates STAT proteins (signal transducers and activators of transcription) and plays a key role in interferon-alpha/beta and interferon-gamma signal transduction. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2016],

Research Area

Image Data



T47D cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with primary antibody 1:1000. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

Note

For research use only.