

Product Name: Rictor Rabbit Monoclonal Antibody
Catalog #: AMRe21154



Summary

Production Name	Rictor Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC,IF,IP,ELISA
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG,Kappa
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
Purification	Protein A

Immunogen

Gene Name	RICTOR
Alternative Names	Rapamycin-insensitive companion of mTOR;AVO3 homolog;hAVO3;
Gene ID	253260.0
SwissProt ID	Q6R327.

Application

Dilution Ratio	IHC 1:200-1:1000;WB 1:2000-1:10000;IF 1:200-1:1000;ELISA 1:5000-1:20000;IP 1:50-1:200;
Molecular Weight	Calculated MW:192kD;Observed MW:192kD

Background

Cell localization:Cytoplasm.RICTOR and MTOR (FRAP1; MIM 601231) are components of a protein complex that integrates


Product Name: Rictor Rabbit Monoclonal Antibody
Catalog #: AMRe21154



nutrient- and growth factor-derived signals to regulate cell growth (Sarbasov et al., 2004 [PubMed 15268862]).[supplied by OMIM, Mar 2008],

Research Area

Image Data



C2C12 whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with primary antibody(1:1000).
The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

Note

For research use only.