

Product Name: FAM160A2 Rabbit Monoclonal Antibody
Catalog #: AMRe21203



Summary

Production Name	FAM160A2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG,Kappa
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
Purification	Protein A

Immunogen

Gene Name	FHIP1B
Alternative Names	FHIP1B;C11orf56;FAM160A2;KIAA1759;FHF complex subunit HOOK-interacting protein 1B ;FHIP1B;FTS- and Hook-interacting protein;FHIP;
Gene ID	84067.0
SwissProt ID	Q8N612.

Application

Dilution Ratio	WB 1:1000-5000
Molecular Weight	Calculated MW;;Observed MW:107kD

Background

The protein encoded by this gene is part of the FTS/Hook/FHIP (FHF) complex, which can interact with members of the


Product Name: FAM160A2 Rabbit Monoclonal Antibody
Catalog #: AMRe21203



homotypic vesicular protein sorting (HOPS) complex. This interaction suggests that the encoded protein is involved in vesicle trafficking. [provided by RefSeq, Dec 2016]

Research Area

Image Data



Western Blot analysis of L02 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-FAM160A2 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

Note

For research use only.