

Product Name: C3orf38 Rabbit Monoclonal Antibody
Catalog #: AMRe21204



Summary

Production Name	C3orf38 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG,Kappa
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
Purification	Protein A

Immunogen

Gene Name	C3orf38
Alternative Names	C3orf38;Uncharacterized protein C3orf38;
Gene ID	285237.0
SwissProt ID	Q5JPI3.

Application

Dilution Ratio	WB 1:1000-5000
Molecular Weight	Calculated MW::Observed MW:36kD

Background


Involved in positive regulation of apoptotic process. Located in nucleus. [provided by Alliance of Genome Resources, Feb 2025]

Product Name: C3orf38 Rabbit Monoclonal Antibody
Catalog #: AMRe21204



Research Area

Image Data



Western Blot analysis of L02 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-C3orf38 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

Note

For research use only.