

---

**Product Name: TAGAP Rabbit Monoclonal Antibody****Catalog #: AMRe21240**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,FC,IP
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG,Kappa
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	0.2mg/ml. The concentration of this product may be batch-dependent.
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%protective protein
<b>Purification</b>	Protein A

**Application**

<b>Dilution Ratio</b>	WB 1:1000-1:5000,FC 1:100-1:300,IP 1:50-1:100
<b>Molecular Weight</b>	Calculated MW;;Observed MW:80kD

**Antigen Information**

<b>Gene Name</b>	TAGAP;TAGAP1;FKSG15
<b>Alternative Names</b>	TAGAP;TAGAP1;FKSG15;T-cell activation Rho GTPase-activating protein ;T-cell activation GTPase-activating protein;
<b>Gene ID</b>	117289.0
<b>SwissProt ID</b>	Q8N103
<b>Immunogen</b>	Recombinant protein of human TAGAP


**Background**

This gene encodes a member of the Rho GTPase-activator protein superfamily. The encoded protein may function as a Rho GTPase-activating protein. Alterations in this gene may be associated with several diseases, including rheumatoid arthritis,

celiac disease, and multiple sclerosis. Alternate splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2013],

## Research Area

## Image Data



Western Blot analysis of HEK-293 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-TAGAP rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.