

Product Name: DNA PKcs Rabbit Monoclonal Antibody
Catalog #: AMRe21286



Summary

Production Name	DNA PKcs Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC,IF,IP,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG,Kappa
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
Purification	Protein A

Immunogen

Gene Name	PRKDC
Alternative Names	PRKDC;HYRC;HYRC1;DNA-dependent protein kinase catalytic subunit;DNA-PK catalytic subunit;DNA-PKcs;DNPK1;p460
Gene ID	5591.0
SwissProt ID	P78527.

Application

Dilution Ratio	IHC 1:200-1:1000;WB 1:2000-1:10000;IF 1:200-1:1000;ELISA 1:5000-1:20000;IP 1:50-1:200;
Molecular Weight	Calculated MW:469kD;Observed MW:469kD

Background


Product Name: DNA PKcs Rabbit Monoclonal Antibody
Catalog #: AMRe21286



Cell localization:Nucleus.This gene encodes the catalytic subunit of the DNA-dependent protein kinase (DNA-PK). It functions with the Ku70/Ku80 heterodimer protein in DNA double strand break repair and recombination. The protein encoded is a member of the PI3/PI4-kinase family.[provided by RefSeq, Jul 2010],

Research Area

Image Data



HT-29 whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with primary antibody(1:1000).
The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

Note

For research use only.