

Product Name: ACE2 Rabbit Monoclonal Antibody
Catalog #: AMRe21429



Summary

Production Name	ACE2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC,IF,IP,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG,Kappa
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
Purification	Protein A

Immunogen

Gene Name	ACE2 Angiotensin-converting enzyme 2;ACE-related carboxypeptidase;Angiotensin-
Alternative Names	converting enzyme homolog;ACEH;Metalloprotease MPROT15;[Cleaved into: Processed angiotensin-converting enzyme 2]
Gene ID	59272.0
SwissProt ID	Q9BYF1.

Application

Dilution Ratio	IHC 1:1000-1:4000;WB 1:1000-1:5000;IF 1:200-1:1000;ELISA 1:5000-1:20000;IP 1:50- 1:200;
Molecular Weight	Calculated MW:130kD;Observed MW:130kD

Product Name: ACE2 Rabbit Monoclonal Antibody
Catalog #: AMRe21429




Background

Cell localization:Membrane.

Research Area

Image Data



Western blot analysis of lysates from CACO-2

cells, using ACE2 Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

Note

For research use only.