

Summary

Production Name	Hsp90 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC,IF,IP,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG,Kappa
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%protective protein
Purification	Protein A

Immunogen

Gene Name	HSP90AB1
Alternative Names	HSP90AB1;HSP90B;HSPC2;HSPCB;Heat shock protein HSP 90-beta;HSP 90;Heat shock 84 kDa;HSP 84;HSP84
Gene ID	3326.0
SwissProt ID	P08238.

Application

Dilution Ratio	IHC 1:200-1000;WB 1:1000-5000;IF 1:200-1000;ELISA 1:5000-20000;IP 1:50-200
Molecular Weight	Calculated MW:85kD;Observed MW:85kD

Background


Product Name: Hsp90 Rabbit Monoclonal Antibody
Catalog #: AMRe21616




Cell localization: Cytoplasm, Nuclear. This gene encodes a member of the heat shock protein 90 family; these proteins are involved in signal transduction, protein folding and degradation and morphological evolution. This gene encodes the constitutive form of the cytosolic 90 kDa heat-shock protein and is thought to play a role in gastric apoptosis and inflammation. Alternative splicing results in multiple transcript variants. Pseudogenes have been identified on multiple chromosomes. [provided by RefSeq, Dec 2012],

Research Area

Image Data



Western blot analysis of lysates from PC12 cells, using Hsp90 Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.



Immunohistochemical analysis of paraffin-embedded Human testis tissue. 1, primary Antibody was diluted at 1:200 (4°C, overnight). 2, EDTA pH 9.0 was used for antibody retrieval (>98°C, 20min). 3, Secondary antibody was diluted at 1:200 (room temperature, 30min).

Note

For research use only.

Product Name: Hsp90 Rabbit Monoclonal Antibody
Catalog #: AMRe21616

