

Product Name: C10orf63 Rabbit Monoclonal Antibody**Catalog #: AMRe83692**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	29 kDa

Antigen Information

Gene Name	C10orf63
Alternative Names	C10orf63; ENKUR; Enkurin;;Enkurin
Gene ID	
SwissProt ID	Q8TC29
Immunogen	A synthesized peptide derived from human Enkurin


Background

Adapter that functions to localize a calcium-sensitive signal transduction machinery in sperm to a calcium-permeable ion channel. Microtubule inner protein (MIP) part of the dynein-decorated doublet microtubules (DMTs) in cilia axoneme, which is

required for motile cilia beating.

Research Area

Image Data



Western blot analysis of C10orf63 expression in Mouse testis lysate.