

---

**Product Name: USP39 Rabbit Monoclonal Antibody****Catalog #: AMRe83696**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,ICC/IF,ICC
<b>Reactivity</b>	Human,Mouse,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
<b>Molecular Weight</b>	65 kDa

**Antigen Information**

<b>Gene Name</b>	USP39
<b>Alternative Names</b>	Inactive ubiquitin specific peptidase 39; SAD1; snRNP ASSEMBLY DEFECTIVE 1; SNRNP65; Ubiquitin specific peptidase 39; USP39;;USP39
<b>Gene ID</b>	
<b>SwissProt ID</b>	Q53GS9
<b>Immunogen</b>	A synthesized peptide derived from human USP39


**Background**

Deubiquitinating enzyme that plays a role in many cellular processes including cellular antiviral response, epithelial

morphogenesis, DNA repair or B-cell development. Plays a role in pre-mRNA splicing as a component of the U4/U6-U5 tri-snRNP, one of the building blocks of the precatalytic spliceosome.

## Research Area

## Image Data



Western blot analysis of USP39 expression in HeLa cell lysate.