
Product Name: Flotillin1 Rabbit Monoclonal Antibody**Catalog #: AMRe83724**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.38mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	Calculated MW: 47 kDa ; Observed MW: 48 kDa

Antigen Information

Gene Name	Flotillin1
Alternative Names	Flotillin1 ; Flotillin 1; Reggie 2; FLOT1; Integral membrane component of caveolae; ;Flotillin 1
Gene ID	
SwissProt ID	O75955
Immunogen	A synthesized peptide derived from human Flotillin 1


Background

May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles.

Research Area

Image Data

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of Flotillin 1 expression in HeLa cell lysate.

