

Summary

Production Name	TCEB2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	TCEB2
Alternative Names	EloB; Elongin B; SIII; SIII p18; Tceb2;;Elongin B
Gene ID	
SwissProt ID	Q15370.A synthesized peptide derived from human Elongin B

Application

Dilution Ratio	WB:1:1000-1:2000
Molecular Weight	Calculated MW: 13 kDa ; Observed MW: 16 kDa

Background


Product Name: TCEB2 Rabbit Monoclonal Antibody
Catalog #: AMRe83729




SIII, also known as elongin, is a general transcription elongation factor that increases the RNA polymerase II transcription elongation past template-encoded arresting sites. Subunit A is transcriptionally active and its transcription activity is strongly enhanced by binding to the dimeric complex of the SIII regulatory subunits B and C (elongin BC complex)

Research Area

Image Data



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of TCEB2 expression in Jurkat cell lysate.

Note

For research use only.