
Product Name: RAIDD Rabbit Monoclonal Antibody**Catalog #: AMRe83731**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ICC,FC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
Molecular Weight	23 kDa

Antigen Information


Gene Name	RAIDD
Alternative Names	CRADD;MGC9163;RAIDD;Death adaptor molecule RAIDD;Death domain containing protein CRADD;;CRADD
Gene ID	
SwissProt ID	P78560
Immunogen	A synthesized peptide derived from human CRADD

Background

Adapter protein that associates with PIDD1 and the caspase CASP2 to form the PIDDosome, a complex that activates CASP2 and triggers apoptosis.

Research Area

Image Data



Western blot analysis of RAIDD expression in HeLa cell lysate.