

Product Name: alpha/beta Synuclein Rabbit Monoclonal Antibody
Catalog #: AMRe83742

Summary

Production Name	alpha/beta Synuclein Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC,FC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	alpha/beta Synuclein
Alternative Names	alpha synuclein; beta synuclein ; NACP; PD1 ; SNCA; SNCB;;alpha/beta Synuclein
Gene ID	
SwissProt ID	P37840/Q16143. A synthesized peptide derived from human alpha Synuclein

Application

Dilution Ratio	WB:1:1000-1:2000,ICC:1:50-1:200,FC:1:20-1:100
Molecular Weight	Calculated MW: 14 kDa ; Observed MW: 18 kDa


Background

Product Name: alpha/beta Synuclein Rabbit Monoclonal Antibody
Catalog #: AMRe83742

Alpha-synuclein is a member of the synuclein family, which also includes beta- and gamma-synuclein. Synucleins are abundantly expressed in the brain and alpha- and beta-synuclein inhibit phospholipase D2 selectively.

Research Area

Image Data



Western blot analysis of alpha + beta Synuclein expression in Mouse brain tissue lysate.

Note

For research use only.