

Product Name: LAP3 Rabbit Monoclonal Antibody
Catalog #: AMRe83762



Summary

Production Name	LAP3 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	LAP3
Alternative Names	LAP; Lap3; LAPEP; Leucine aminopeptidase 3; PEPS; Peptidase S;;LAP3
Gene ID	
SwissProt ID	P28838.A synthesized peptide derived from human LAP3

Application

Dilution Ratio	WB:1:1000-1:2000,IHC:1:100-1:200
Molecular Weight	Calculated MW: 56 kDa ; Observed MW: 52 kDa

Background


Product Name: LAP3 Rabbit Monoclonal Antibody
Catalog #: AMRe83762




Presumably involved in the processing and regular turnover of intracellular proteins. Catalyzes the removal of unsubstituted N-terminal amino acids from various peptides.

Research Area

Image Data



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of LAP3 expression in 293T cell lysate.

Note

For research use only.