

---

**Product Name: Tripeptidyl peptidase 2 Rabbit Monoclonal Antibody****Catalog #: AMRe83763**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,ICC/IF,ICC,FC
<b>Reactivity</b>	Human,Mouse,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	0.55mg/ml. The concentration of this product may be batch-dependent.
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
<b>Molecular Weight</b>	138 kDa

**Antigen Information**

<b>Gene Name</b>	Tripeptidyl peptidase 2
<b>Alternative Names</b>	TPP II; TPP2; Tripeptidyl aminopeptidase; Tripeptidyl peptidase 2; Tripeptidyl peptidase II;;TPP2
<b>Gene ID</b>	
<b>SwissProt ID</b>	P29144
<b>Immunogen</b>	A synthesized peptide derived from human TPP2


**Background**

Tripeptidyl peptidase II is a component of the proteolytic cascade acting downstream of the 26S proteasome in the ubiquitin-


proteasome pathway. It may be able to complement the 26S proteasome function to some extent under conditions in which the latter is inhibited.

## Research Area


## Image Data



Western blot analysis of Tripeptidyl peptidase II expression in HepG2 cell lysate.



Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunofluorescent analysis using the Antibody at 1:50 dilution.