

Product Name: RALBP1 Rabbit Monoclonal Antibody
Catalog #: AMRe83807



Summary

Production Name	RALBP1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC,FC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	RALBP1
Alternative Names	RALBP1; RIP1; RLIP1; RLIP76;;RalBP1
Gene ID	
SwissProt ID	Q15311.A synthesized peptide derived from human RalBP1

Application

Dilution Ratio	WB:1:1000-1:5000,IHC:1:100-1:200,FC:1:20-1:100
Molecular Weight	Calculated MW: 76 kDa ; Observed MW: 95 kDa

Background


Product Name: RALBP1 Rabbit Monoclonal Antibody
Catalog #: AMRe83807




Multifunctional protein that functions as a downstream effector of RALA and RALB. As a GTPase-activating protein/GAP can inactivate CDC42 and RAC1 by stimulating their GTPase activity. As part of the Ral signaling pathway, may also regulate ligand-dependent EGF and insulin receptors-mediated endocytosis.

Research Area

Image Data



Western blot analysis of RALBP1 expression in HeLa cell lysate.



Immunohistochemical analysis of paraffin-embedded human breast, using RALBP1 Antibody.

Note

For research use only.