
Product Name: AKAP14 Rabbit Monoclonal Antibody**Catalog #: AMRe83817**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,ICC 1:50-1:200,IP 1:20-1:50
Molecular Weight	23 kDa

Antigen Information


Gene Name	AKAP14
Alternative Names	A-kinase anchor protein 14; A-kinase anchor protein 14; A-kinase anchor protein 28 kDa; AKAP14; AKAP28;;AKAP 14
Gene ID	
SwissProt ID	Q86UN6
Immunogen	A synthesized peptide derived from human AKAP 14

Background

Binds to type II regulatory subunits of protein kinase A and anchors/targets them.

Research Area

Image Data



Western blot analysis of AKAP14 expression in (1) Jurkat cell lysate; (2) A549 cell lysate.