Product Name: BRG1 Rabbit Monoclonal Antibody

Catalog #: AMRe83832



Summary

Production Name BRG1 Rabbit Monoclonal Antibody

Description Rabbit Monoclonal antibody

Host Rabbit

Application WB,IP,IHC,ICC,IF **Reactivity** Human,Mouse,Rat

Performance

ConjugationUnconjugatedModificationUnmodified

Isotype IgG

Clonality Monoclonal Form Liquid

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw Storage

cycles.

Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% **Buffer**

glycerol.

Purification Affinity Purification

Immunogen

Gene Name BRG1

ATP-dependent helicase SMARCA4; BRG1-associated factor 190A; BAF190A; Mitotic

Alternative Names growth and transcription activator; Protein BRG-1; Protein brahma homolog 1; SNF2-

beta;;SMARCA4

Gene ID

SwissProt ID P51532.A synthesized peptide derived from human SMARCA4

Application

Dilution Ratio WB:1:1000-1:5000,IP:ChIP:1:20-1:50,IHC:1:100-1:200;IF:1:50-1:200,ICC:1:50-1:200

Molecular Weight 185 kDa

Product Name: BRG1 Rabbit Monoclonal Antibody

Catalog #: AMRe83832




Background


Involved in transcriptional activation and repression of select genes by chromatin remodeling (alteration of DNA-nucleosome topology). Component of SWI/SNF chromatin remodeling complexes that carry out key enzymatic activities, changing chromatin structure by altering DNA-histone contacts within a nucleosome in an ATP-dependent manner.

Research Area

Image Data




All lanes use the Antibody at 1:5K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:5K dilution for 1 hour at room temperature.

Product Name: BRG1 Rabbit Monoclonal Antibody Catalog #: AMRe83832

C EnkiLife



Western blot analysis of BRG1 expression in HeLa cell lysate.

Note

For research use only.