

Summary

Production Name	IRF6 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC,IF
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	IRF6
Alternative Names	Interferon regulatory factor 6; IRF 6; LPS;O6; PIT; PPS; PPS1; VWS; VWS1;;IRF6
Gene ID	
SwissProt ID	O14896.A synthesized peptide derived from human IRF6

Application

Dilution Ratio	WB:1:1000-1:2000,IHC:1:100-1:200;IF:1:50-1:200,ICC:1:50-1:200
Molecular Weight	53 kDa

Background


Product Name: IRF6 Rabbit Monoclonal Antibody
Catalog #: AMRe83836




Probable DNA-binding transcriptional activator. Key determinant of the keratinocyte proliferation-differentiation switch involved in appropriate epidermal development.

Research Area


Image Data



Western blot analysis of IRF6 expression in Jurkat cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Product Name: IRF6 Rabbit Monoclonal Antibody
Catalog #: AMRe83836



Note

For research use only.