
Product Name: Hsp20 Rabbit Monoclonal Antibody**Catalog #: AMRe83870**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,FC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:5000,IHC 1:100-1:200,FC 1:20-1:100
Molecular Weight	17 kDa

Antigen Information


Gene Name	Hsp20
Alternative Names	HSPB6; HEL55; Hsp20; Heat-shock 20 kDa like-protein p20;;Hsp20
Gene ID	
SwissProt ID	O14558
Immunogen	A synthesized peptide derived from human Hsp20

Background

Small heat shock protein which functions as a molecular chaperone probably maintaining denatured proteins in a folding-competent state. Seems to have versatile functions in various biological processes.

Research Area

Image Data



Western blot analysis of Hsp20 expression in mouse skeletal muscle lysate.