

Product Name: GC/VDBP Rabbit Monoclonal Antibody**Catalog #: AMRe83906**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ICC,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
Molecular Weight	53 kDa

Antigen Information


Gene Name	GC/VDBP
Alternative Names	DBP; DBP/GC; GRD3; VDBG; VDBP;;Vitamin D binding protein
Gene ID	
SwissProt ID	P02774
Immunogen	A synthesized peptide derived from human Vitamin D binding protein

Background

Involved in vitamin D transport and storage, scavenging of extracellular G-actin, enhancement of the chemotactic activity of C5 alpha for neutrophils in inflammation and macrophage activation.

Research Area

Image Data



Western blot analysis of GC/VDBP expression in Human plasma lysate.