

Product Name: Cyclophilin A Rabbit Monoclonal Antibody**Catalog #: AMRe83911**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200
Molecular Weight	18 kDa

Antigen Information

Gene Name	Cyclophilin A
Alternative Names	CyclophilinA; Cyclosporin A binding protein; CYPA; CYPH; Peptidyl-prolyl cis-trans isomerase A; Rotamase A; PPlase A; PPIA; Peptidyl-prolyl cis-trans isomerase A;;PPIA
Gene ID	
SwissProt ID	P62937
Immunogen	A synthesized peptide derived from human PPIA


Background

PPlases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in

oligopeptides.

Research Area

Image Data



Western blot analysis of Cyclophilin A expression in HepG2 cell lysate.