

Product Name: Endothelin B Receptor Rabbit Monoclonal Antibody**Catalog #: AMRe83912**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000, IP 1:20-1:50
Molecular Weight	50 kDa

Antigen Information

Gene Name	Endothelin B Receptor
Alternative Names	ABCDS; Ednra; EDNRB; Endothelin B receptor precursor; ETB; ETBR; ETRB; Hirschsprung disease;;EDNRB
Gene ID	
SwissProt ID	P24530
Immunogen	A synthesized peptide derived from human EDNRB


Background

Non-specific receptor for endothelin 1, 2, and 3. Mediates its action by association with G proteins that activate a

phosphatidylinositol-calcium second messenger system.

Research Area

Image Data



Western blot analysis of Endothelin B Receptor expression in 293T cell lysate.