

---

**Product Name: SFRP1 Rabbit Monoclonal Antibody****Catalog #: AMRe83935**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,ICC/IF,ICC
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	0.38mg/ml. The concentration of this product may be batch-dependent.
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
<b>Molecular Weight</b>	35 kDa

**Antigen Information**


<b>Gene Name</b>	SFRP1 Frizzled related protein 1; FRP1; FrzA; SARP2; Secreted apoptosis related protein 2; Secreted
<b>Alternative Names</b>	frizzled related protein 1; SFRP1; ;FRP 1
<b>Gene ID</b>	
<b>SwissProt ID</b>	Q8N474
<b>Immunogen</b>	A synthesized peptide derived from human FRP 1

**Background**


Soluble frizzled-related proteins (sFRPS) function as modulators of Wnt signaling through direct interaction with Wnts. They have a role in regulating cell growth and differentiation in specific cell types. SFRP1 decreases intracellular beta-catenin levels.

## Research Area

## Image Data



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of SFRP1 expression in A375 cell lysate.