
Product Name: GCET1 Rabbit Monoclonal Antibody**Catalog #: AMRe83943**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	47 kDa

Antigen Information


Gene Name	GCET1
Alternative Names	Centerin; GCET1; Seprin A9; SERPINA11; SERPINA11b; Serpina9;;Serpina A9
Gene ID	
SwissProt ID	Q86WD7
Immunogen	A synthesized peptide derived from human Serpin A9

Background

Protease inhibitor that inhibits trypsin and trypsin-like serine proteases (in vitro). Inhibits plasmin and thrombin with lower efficiency (in vitro).

Research Area

Image Data



Western blot analysis of GCET1 expression in MCF7 cell lysate.