
Product Name: Reelin Rabbit Monoclonal Antibody**Catalog #: AMRe83986**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	388 kDa

Antigen Information

Gene Name	Reelin
Alternative Names	LIS2; PRO1598; Reeler; Reelin; RELN; RL;;Reelin
Gene ID	
SwissProt ID	P78509
Immunogen	A synthesized peptide derived from human Reelin


Background

Extracellular matrix serine protease that plays a role in layering of neurons in the cerebral cortex and cerebellum. Regulates microtubule function in neurons and neuronal migration. Affects migration of sympathetic preganglionic neurons in the spinal

cord, where it seems to act as a barrier to neuronal migration.

Research Area

Image Data



Western blot analysis of Reelin expression in HepG2 cell lysate.