
Product Name: ProteinC Rabbit Monoclonal Antibody**Catalog #: AMRe83988**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	52 kDa

Antigen Information

Gene Name	ProteinC
Alternative Names	Activation peptide; Anticoagulant protein C; APC; Autoprothrombin IIA; Blood coagulation factor XIV; PROC1; THPH3; THPH4; Vitamin K dependent protein C; ;PROC
Gene ID	
SwissProt ID	P04070
Immunogen	A synthesized peptide derived from human PROC


Background

Protein C is a vitamin K-dependent serine protease that regulates blood coagulation by inactivating factors Va and VIIIa in the

presence of calcium ions and phospholipids. Exerts a protective effect on the endothelial cell barrier function.

Research Area

Image Data



Western blot analysis of Protein C expression in HepG2 cell lysate.