

Product Name: Chromogranin B Rabbit Monoclonal Antibody**Catalog #: AMRe83989**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	Calculated MW: 78 kDa ; Observed MW: 125 kDa

Antigen Information

Gene Name	Chromogranin B
Alternative Names	CCB peptide; CgB; Chgb; ChromograninB; SCG1; Secretogranin 1; Secretogranin B; Sgl;;Secretogranin 1
Gene ID	
SwissProt ID	P05060
Immunogen	A synthesized peptide derived from human Secretogranin 1


Background

Secretogranin-1 is a neuroendocrine secretory granule protein, which may be the precursor for other biologically active

peptides.

Research Area

Image Data



Western blot analysis of Chromogranin B expression in (1) HeLa cell lysate; (2) Raw 264.7 cell lysate.