
Product Name: Dynein intermediate chain 1 Rabbit Monoclonal Antibody**Catalog #: AMRe83991**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	79 kDa

Antigen Information


Gene Name	Dynein intermediate chain 1
Alternative Names	DNAI1; ICS; Immotile cilia syndrome 1; ICS1; CILD1; PCD;;DNAI1
Gene ID	
SwissProt ID	Q9UI46
Immunogen	A synthesized peptide derived from human DNAI1

Background

Part of the dynein complex of respiratory cilia.

Research Area

Image Data



Western blot analysis of Dynein intermediate chain 1 expression in mouse testis lysate.