

Product Name: PI3 Kinase p110 delta Rabbit Monoclonal Antibody
Catalog #: AMRe84012

Summary

Production Name	PI3 Kinase p110 delta Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	PI3 Kinase p110 delta
Alternative Names	APDS; GRB1; IMD14; p110dp85a; p85-ALPHA; Phosphoinositide 3 kinase B; Phosphoinositide 3 kinase C; p37delta; PI3 kinase p110 subunit delta; PI3Kdelta; Pik3cd; PIK3R1; PK3CD; PtdIns 3 kinase p110;;PI3 kinase p110 delta
Gene ID	
SwissProt ID	O00329.A synthesized peptide derived from human PI3 kinase p110 delta

Application

Dilution Ratio	WB:1:1000-1:2000,ICC:1:50-1:200
Molecular Weight	119 kDa


Product Name: PI3 Kinase p110 delta Rabbit Monoclonal Antibody
Catalog #: AMRe84012

Background


Phosphoinositide-3-kinase (PI3K) phosphorylates phosphatidylinositol (PI) and its phosphorylated derivatives at position 3 of the inositol ring to produce 3-phosphoinositides.

Research Area

Image Data




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.


Product Name: PI3 Kinase p110 delta Rabbit Monoclonal Antibody
Catalog #: AMRe84012



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.




Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunofluorescent analysis using the Antibody at 1:50 dilution.

Product Name: PI3 Kinase p110 delta Rabbit Monoclonal
Antibody
Catalog #: AMRe84012



Immunofluorescent analysis using the Antibody at 1:50 dilution.

Note

For research use only.