
Product Name: RANTES(CCL5) Rabbit Monoclonal Antibody**Catalog #: AMRe84015**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,ICC 1:50-1:200
Molecular Weight	10 kDa

Antigen Information


Gene Name	RANTES(CCL5)
Alternative Names	SISd; eoCP; SCYA5; RANTES; TCP228; D17S136E; SIS-delta;;CCL5
Gene ID	
SwissProt ID	P13501
Immunogen	A synthesized peptide derived from human CCL5

Background

Chemoattractant for blood monocytes, memory T-helper cells and eosinophils. Causes the release of histamine from basophils and activates eosinophils. May activate several chemokine receptors including CCR1, CCR3, CCR4 and CCR5.

Research Area

Image Data



Western blot analysis of RANTES expression in RANTES recombinant protein lysate.