

Product Name: PDE4D Rabbit Monoclonal Antibody**Catalog #: AMRe84026**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	91 kDa

Antigen Information


Gene Name	PDE4D
Alternative Names	DPDE3; HSPDE4D; PDE43; Pde4d; PDE4DN2; STRK1;;PDE4D
Gene ID	
SwissProt ID	Q08499
Immunogen	A synthesized peptide derived from human PDE4D

Background

Hydrolyzes the second messenger cAMP, which is a key regulator of many important physiological processes.

Research Area

Image Data



Western blot analysis of PDE4D expression in K562 cell lysate.