

**Product Name: MST4 Rabbit Monoclonal Antibody****Catalog #: AMRe84037**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC
<b>Reactivity</b>	Human,Mouse,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:1000-1:2000,IHC 1:100-1:200
<b>Molecular Weight</b>	Calculated MW: 47 kDa ; Observed MW: 50 kDa

**Antigen Information**


<b>Gene Name</b>	MST4
<b>Alternative Names</b>	MASK; Mst4; STE20 like kinase MST4; STK26;;MST4
<b>Gene ID</b>	
<b>SwissProt ID</b>	Q9P289
<b>Immunogen</b>	A synthesized peptide derived from human MST4

**Background**

Serine/threonine-protein kinase that acts as a mediator of cell growth.Modulates apoptosis.In association with STK24 negatively regulates Golgi reorientation in polarized cell migration upon RHO activation.

## Research Area

## Image Data



Western blot analysis of MST4 expression in HeLa cell lysate.