

Product Name: MS4A14 Rabbit Monoclonal Antibody
Catalog #: AMRe84038



Summary

Production Name	MS4A14 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	MS4A14
Alternative Names	MS4A13; MS4A16; NYDSP21;;MS4A14
Gene ID	
SwissProt ID	Q96JA4. A synthesized peptide derived from human MS4A14

Application

Dilution Ratio	WB:1:1000-1:2000
Molecular Weight	Calculated MW: 77 kDa ; Observed MW: 80 kDa

Background


May be involved in signal transduction as a component of a multimeric receptor complex.

Product Name: MS4A14 Rabbit Monoclonal Antibody
Catalog #: AMRe84038



Research Area

Image Data



Western blot analysis of MS4A14 expression in MCF7 cell lysate.

Note

For research use only.