
Product Name: GALT Rabbit Monoclonal Antibody**Catalog #: AMRe84039**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,ICC 1:50-1:200
Molecular Weight	43 kDa

Antigen Information


Gene Name	GALT
Alternative Names	GALT;;GALT
Gene ID	
SwissProt ID	P07902
Immunogen	A synthesized peptide derived from human GALT

Background

Plays an important role in galactose metabolism.

Research Area

Image Data



Western blot analysis of GALT expression in K562 cell lysate.