
Product Name: KCTD9 Rabbit Monoclonal Antibody**Catalog #: AMRe84040**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,ICC 1:50-1:200
Molecular Weight	43 kDa

Antigen Information


Gene Name	KCTD9
Alternative Names	BTBD27;;KCTD9
Gene ID	
SwissProt ID	Q7L273
Immunogen	A synthesized peptide derived from human KCTD9

Background

Substrate-specific adapter of a BCR (BTB-CUL3-RBX1) E3 ubiquitin-protein ligase complex, which mediates the ubiquitination of target proteins, leading to their degradation by the proteasome.

Research Area

Image Data



Western blot analysis of KCTD9 expression in HepG2 cell lysate.