

Summary

Production Name	AAMP Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC,FC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	AAMP
Alternative Names	AAMP;;AAMP
Gene ID	
SwissProt ID	Q13685.A synthesized peptide derived from human AAMP

Application

Dilution Ratio	WB:1:1000-1:2000,ICC:1:50-1:200,FC:1:20-1:100
Molecular Weight	Calculated MW: 47 kDa ; Observed MW: 52 kDa

Background


Product Name: AAMP Rabbit Monoclonal Antibody
Catalog #: AMRe84042




Plays a role in angiogenesis and cell migration. In smooth muscle cell migration, may act through the RhoA pathway.

Research Area


Image Data



Western blot analysis of AAMP expression in (1) A375 lysate; (2) MCF7 cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

Product Name: AAMP Rabbit Monoclonal Antibody
Catalog #: AMRe84042



Note

For research use only.