

Product Name: DAZAP1 Rabbit Monoclonal Antibody
Catalog #: AMRe84044



Summary

Production Name	DAZAP1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	DAZAP1
Alternative Names	DAZ associated protein 1; Dazap1;;DAZ associated protein 1
Gene ID	
SwissProt ID	Q96EP5. A synthesized peptide derived from human DAZ associated protein 1

Application

Dilution Ratio	WB:1:1000-1:2000
Molecular Weight	43 kDa

Background


RNA-binding protein, which may be required during spermatogenesis.

Product Name: DAZAP1 Rabbit Monoclonal Antibody
Catalog #: AMRe84044



Research Area

Image Data



Western blot analysis of DAZAP1 expression in Jurkat cell lysate.

Note

For research use only.