

Summary

Production Name	CECR5 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	CECR5
Alternative Names	Cecr5;;HDHD5
Gene ID	
SwissProt ID	Q9BXW7.A synthesized peptide derived from human HDHD5

Application

Dilution Ratio	WB:1:1000-1:2000, ICC:1:50-1:200
Molecular Weight	46 kDa

Background


Product Name: CECR5 Rabbit Monoclonal Antibody
Catalog #: AMRe84045




Candidate gene for the Cat Eye Syndrome (CES), a developmental disorder associated with the duplication of a 2 Mb region of 22q11.2. Duplication usually takes in the form of a supernumerary bisatellited isodicentric chromosome, resulting in four copies of the region (represents an inv dup22(q11)). CES is characterized clinically by the combination of coloboma of the iris and anal atresia with fistula, downslanting palpebral fissures, preauricular tags and/or pits, frequent occurrence of heart and renal malformations, and normal or near-normal mental development.

Research Area

Image Data



Western blot analysis of CECR5 expression in HepG2 cell lysate.



All lanes use the Antibody at 1:1k dilution for 1 hour at room temperature.

Note

For research use only.