

---

**Product Name: TMEM192 Rabbit Monoclonal Antibody****Catalog #: AMRe84053**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,ICC,FC
<b>Reactivity</b>	Human,Mouse,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	0.59mg/ml. The concentration of this product may be batch-dependent.
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:1000-1:2000,ICC 1:50-1:200,FC 1:20-1:100
<b>Molecular Weight</b>	Calculated MW: 31 kDa ; Observed MW: 29 kDa

**Antigen Information**


<b>Gene Name</b>	TMEM192
<b>Alternative Names</b>	TMEM192; Transmembrane protein 192;;TMEM192
<b>Gene ID</b>	
<b>SwissProt ID</b>	Q8IY95
<b>Immunogen</b>	A synthesized peptide derived from human TMEM192

**Background**


Enables protein homodimerization activity. Located in several cellular components, including late endosome; lysosomal membrane; and perinuclear region of cytoplasm.

## Research Area


## Image Data



Western blot analysis of TMEM192 expression in U87-MG cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.