
Product Name: ACAP2 Rabbit Monoclonal Antibody**Catalog #: AMRe84063**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.59mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	88 kDa

Antigen Information


Gene Name	ACAP2
Alternative Names	ACAP2; Centaurin beta2; CENTB2; Cntb2;;ACAP2
Gene ID	
SwissProt ID	Q15057
Immunogen	A synthesized peptide derived from human ACAP2

Background


GTPase-activating protein (GAP) for ADP ribosylation factor 6 (ARF6).

Research Area


Image Data



Western blot analysis of ACAP2 expression in 293 cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.