

Product Name: PGM1 Rabbit Monoclonal Antibody**Catalog #: AMRe84066**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,ICC 1:50-1:200,IP 1:20-1:50
Molecular Weight	61 kDa

Antigen Information


Gene Name	PGM1
Alternative Names	CDG1T; Glucose phosphomutase 1; GSD14; PGM1; Phosphoglucomutase 1;;PGM 1
Gene ID	
SwissProt ID	P36871
Immunogen	A synthesized peptide derived from human PGM 1

Background

This enzyme participates in both the breakdown and synthesis of glucose.

Research Area

Image Data



Western blot analysis of PGM1 expression in (1) 293 cell lysate; (2) RAW 264.7 cell lysate.