
Product Name: RNF115 Rabbit Monoclonal Antibody**Catalog #: AMRe84068**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	Calculated MW: 34 kDa ; Observed MW: 37 kDa

Antigen Information


Gene Name	RNF115
Alternative Names	BCA2; RNF115;;RNF115
Gene ID	
SwissProt ID	Q9Y4L5
Immunogen	A synthesized peptide derived from human RNF115

Background

E3 ubiquitin-protein ligase that catalyzes the 'Lys-48'- and/or 'Lys-63'-linked polyubiquitination of various substrates and thereby plays a role in a number of signaling pathways including autophagy, innate immunity, cell proliferation and cell death.

Research Area

Image Data



Western blot analysis of RNF115 expression in PC-3 cell lysate.