
Product Name: TAGAP Rabbit Monoclonal Antibody**Catalog #: AMRe84077**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IP 1:20-1:50
Molecular Weight	Calculated MW: 81 kDa ; Observed MW: 65,90 kDa

Antigen Information


Gene Name	TAGAP
Alternative Names	ARHGAP47; IDDM21; Tagap; TAGAP1;;TAGAP
Gene ID	
SwissProt ID	Q8N103
Immunogen	A synthesized peptide derived from human TAGAP

Background

May function as a GTPase-activating protein and may play important roles during T-cell activation.

Research Area

Image Data



Western blot analysis of TAGAP expression in Jurkat cell lysate.